

MS29513-020
O-RING INSTALLATION
1 Each
C/D 1Q89
AVI
DAYTON, OH

RR-56S
RING RETAINING
1 Each
AVI
DAYTON, OH

MS29513-010 ORING 1 EA.
PO. 74475 LN. 1
MFG PARKER NITRILE
LOT 326630 DCN 248396
MIL SPEC MIL-P-5315
WENCOR WEST INC.
CURE 3Q98

971C12K KIT (ORING REPLACEMENT)
1 4 x 4 ZIP LOCK BAG
1 TAG PROCEDURE
1 RR56S RETAINING RING
1 MS29513-020 ORING
INSTALLATION C/D 1Q89
1 010-7657 ORING (POPPET)
C/D 1Q88
AVI
DAYTON, OH

O-RING REPLACEMENT KIT P/N 971C12 INSTRUCTIONS

1. Operate valve to locked open prevent spring loaded component flying apart during disassembly
2. Remove and discard spirolax poppet (D) spring (E), retaining O-ring (C) from valve body
3. Remove poppet assembly, discarding (C) from valve body
4. Remove (C) from valve body